

SAMPLE SUBMISSION GUIDELINES

All plates and tubes should be clearly and unmistakably marked

Concentration

Sequencing primer	5 pmol/ μ l.
Purified templates	Plasmid DNA : 50 ng/ μ l PCR products < 1 kb : 20 ng/ μ l PCR products > 1 kb : 100 ng/ μ l

Volume

For both primer and template deliver at least 10 μ l per sequencing reaction that needs to be performed with it.



VIB Genetic Service Facility
University of Antwerp (Campus Drie Eiken, room V0.07)
Universiteitsplein 1
B-2610 Wilrijk
Email gsf@ua.ac.be
Phone +32 3 265 10 99
Fax +32 3 265 11 00